



## Structure & Ordering Information

Product	Catalog No.	Diameter	Usable Length	Spring Coil Length	Coating Length	Radiopaque Length	Tip Shape	Clip Color
ASAHI Gladius MG	AP14R023S	0.36mm (0.014 inch)	190cm	8.5cm	SLIP-COAT® 41cm	3cm	Straight	
	AP14R323S		300cm					Pre-shape
	AP14R023P		190cm				Pre-shape	
	AP14R323P		300cm					
ASAHI Gladius MG ES	AP14R024S	0.36mm (0.014 inch)	190cm	3.0cm	SLIP-COAT® 10cm	3cm	Straight	
	AP14R324S		300cm					Pre-shape
	AP14R024P		190cm				Pre-shape	
	AP14R324P		300cm					

### ASAHI Gladius MG



### ASAHI Gladius MG ES



ASAHI brand multi wire coil provides torque force, torque response, durability and flexibility.

Not all the shapes and sizes are available in all markets. Please check availability with your sales representative.

Some materials/graphics were designed & produced by Vascular Perspective Ltd. [www.vascularperspectives.com](http://www.vascularperspectives.com).

© 2022 ASAHI INTECC CO., LTD. "ASAHI" and "Gladius" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries. Coated with SLIP-COAT® coating

Manufactured by:

Contact:

**ASAHI INTECC**

**ASAHI INTECC CO.,LTD.**

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan  
Contact phone number: +81-561-86-9101

[www.asahi-intecc.com](http://www.asahi-intecc.com)



for Android



for iOS


ASAHI product information

P22021\_B\_NOCE\_2.0



# ASAHI Gladius MG

PCI GUIDE WIRE



YOUR DREAMS. WOVEN TOGETHER.

# ASAHI Gladius MG

Broaden your approach to CTO cases



## Characteristics



- Narrow loop and durability delivered through the modified distal core design
- High lubricity for smooth manipulation
- High torque response provided by ACT ONE technology.




## Excellent Durability

ASAHI Gladius MG maintains shape and guide wire performance even after crossing the lesion.



Guide wire prolapse examination after lesion crossing simulation test



ASAHI Gladius MG ES




Competitor A



### Torque Response Test


ASAHI Gladius MG with ACT ONE provides quick initial torque and superior torque response without whip motion.

— ASAHI Gladius MG  
— ASAHI Gladius MG ES  
— Competitor A  
- - - Ideal line





### Lubricity Test

ASAHI Gladius MG has excellent lubricity that allows for smooth guide wire manipulation.



### Flexibility Test

ASAHI Gladius MG has a modified distal core balance that minimizes a prolapse occurring in the tip.



— ASAHI Gladius EX  
— ASAHI Gladius MG  
— ASAHI Gladius MG ES  
— Competitor A

\*All the data was obtained by company standardized test, which may differ from industry standardized tests.  
\*All the data does not prove that all devices have exactly the same performance with the samples used for these tests.